

KNEE-TUNNEL (CAMP COVENTRY)

ROUTINE: Camp Coventry

DISTANCE: 40"

BUCKY: Yes

FOCAL SPOT: Small

IMAGING PLATE: 10x12 CR cassette or DR detector

CENTRAL RAY:

1. 40° caudal angle entering the popliteal region and exiting $\frac{1}{2}$ " inferior to the patellar apex.

HELPFUL POSITIONING:

Patient is prone with affected knee flexed 40° with support under foot.

Reviewed 10/01/19

